

F · S · T

Surgical Narrow Forceps - Slim1x2 Teeth/Straight/14.5cm



General Info

Product Number	11023-14
Category	Forceps & Tweezers > Toothed
EAN	4251909501495



Find 11023-14 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	25 g	Net Weight	21 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	14.5 cm		
Tipdimensions	2 x 2.2 mm		
Tips	1x2 Teeth / Narrow		
Tipshape	Straight		