



Surgical Narrow Forceps - Slim1x2 Teeth/Straight/15.5cm



Allgemein

Produktnummer	11023-15
Kategorie	Forceps & Tweezers > Toothed
EAN	4061706024024



Finde 11023-15 auf finescience.com

Technische Informationen

Versandgewicht	30 g	Nettogewicht	28 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	15.5 cm		
Tipdimensions	2 x 2.2 mm		
Tips	1x2 Teeth / Narrow		
Tipshape	Straight		