



Semken Surgical Forceps . 1x2 Teeth/Straight/12.5cm



The Semken forceps are characterized by their slim shape and fine tips. Their anatomical (serrationed) or surgical (1x2 teeth) tips make them ideal for grasping tissue. The serrated grip surface ensures a secure hold during use.

General Info

Product Number	11024-12
Category	Forceps & Tweezers > Toothed
EAN	4251909511739



Find 11024-12 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	14.5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	12.5 cm		
Tipdimensions	1.3 x 2 mm		
Tipwidth	1.3 mm		
Tips	1x2 Teeth		
Tipshape	Straight		