



Russian Forceps - Cupped Teeth/Straight/15cm



These forceps have spoon-shaped tips with serrated edges.

General Info

Product Number	11026-15
Category	Forceps & Tweezers > Toothed
EAN	4251909501525



Find 11026-15 on finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	31 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	15 cm		
Tipdimensions	6 x 11 mm		
Tips	Cupped teeth		
Tipshape	Straight		