



Adson Forceps - 1x2 Teeth/Straight/12cm



General Info

Product Number	11027-12
Category	Forceps & Tweezers > Toothed Most Popular > Forceps & Tweezers
EAN	4251909501532



Find 11027-12 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	20 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	12 cm		
Tipdimensions	1.5 x 2 mm		
Tipthickness	2 mm		
Tipwidth	1.5 mm		
Tips	1x2 Teeth		
Tipshape	Straight		