



Gillies Forceps - 1x2 Teeth/Straight/15cm



General Info

Product Number	11028-15
Category	Forceps & Tweezers > Toothed
EAN	4251909501549



Find 11028-15 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	30 g	Net Weight	22 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	15 cm		
Tipdimensions	1.5 x 2 mm		
Tips	1x2 Teeth & cross serrated		
Tipshape	Straight		