



Surgical Narrow Forceps - 1x2 Teeth/Curved/14cm



Allgemein

Produktnummer	11029-14
Kategorie	Forceps & Tweezers > Toothed
EAN	4251909501556



Finde 11029-14 auf finescience.com

Technische Informationen

Versandgewicht	25 g	Nettogewicht	22 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	14 cm		
Tipdimensions	3 x 2 mm		
Tips	1x2 Teeth / Narrow		
Tipshape	Curved		