



Surgical Narrow Forceps - 1x2 Teeth/Curved/14cm



General Info

Product Number	11029-14
Category	Forceps & Tweezers > Toothed
EAN	4251909501556



Find 11029-14 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	22 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	14 cm		
Tipdimensions	3 x 2 mm		
Tips	1x2 Teeth / Narrow		
Tipshape	Curved		