



Transport Forceps - with Screw/Serrated/Straight/14.5cm



Adjustable screw allows for a set holding pressure.

General Info

Product Number	11031-15
Category	Forceps & Tweezers > Serrated
EAN	4251909501570



Find 11031-15 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	28 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	14.5 cm		
Tipdimensions	8 x 11 mm		
Tips	Serrated		
Tipshape	Straight		