



Foerster Forceps - Serrated/Straight/9.5cm



The octagon grip on these forceps provides a secure hold.

General Info

Product Number	11040-08
Category	Forceps & Tweezers > Serrated
EAN	4251909501594



Find 11040-08 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	10 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	9.5 cm		
Tipdimensions	0.8 x 0.5 mm		
Tipwidth	0.8 mm		
Tips	Serrated		
Tipshape	Straight		