



Foerster Forceps - Serrated/Curved/9.5cm



The octagon grip on these forceps provides a secure hold.

General Info

Product Number	11041-08
Category	Forceps & Tweezers > Serrated
EAN	4251909501600



Find 11041-08 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	9 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	9.5 cm		
Tipdimensions	0.8 x 0.5 mm		
Tipwidth	0.8 mm		
Tips	Serrated		
Tipshape	Curved		