



Foerster Forceps - 1x2 Teeth/Curved/9.5cm



The octagon grip on these forceps provides a secure hold.

General Info

Product Number	11043-08
Category	Forceps & Tweezers > Toothed
EAN	4251909501624



Find 11043-08 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	9 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	9.5 cm		
Tipdimensions	0.8 x 0.5 mm		
Tipwidth	0.8 mm		
Tips	1x2 Teeth		
Tipshape	Curved		