



Graefe Forceps - Serrated/Angled/10cm



Allgemein

Produktnummer	11049-10
Kategorie	Forceps & Tweezers > Serrated Most Popular > Forceps & Tweezers
EAN	4251909501631



Finde 11049-10 auf finescience.com

Technische Informationen

Versandgewicht	10 g	Nettogewicht	7 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	10 cm		
Tipdimensions	0.8 x 0.7 mm		
Tipwidth	0.8 mm		
Tips	Serrated		
Tipshape	Angled		