



Graefe forceps are fine instruments which are available in 0.5 or 0.8 mm tip diameters. The tips are either anatomical (serrated) or surgical (1x2 teeth). Thanks to their wide range of variants, they have established themselves in many different applications.

General Info

Product Number	11050-10
Category	Forceps & Tweezers > Serrated Most Popular > Forceps & Tweezers
EAN	4251909501648



Find 11050-10 on
finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	7 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	10 cm		
Tipdimensions	0.8 x 0.7 mm		
Tipwidth	0.8 mm		
Tips	Serrated		
Tipshape	Straight		