



Graefe Forceps - Serrated/Slightly Curved/10cm



General Info

Product Number	11051-10
Category	Forceps & Tweezers > Serrated Most Popular > Forceps & Tweezers
EAN	4061706098551



Find 11051-10 on
finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	7 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	10 cm		
Tipdimensions	0.8 x 0.7 mm		
Tipwidth	0.8 mm		
Tips	Serrated		
Tipshape	Slightly curved		