



Graefe Forceps - 1x2 Teeth/Curved/10cm



Allgemein

Produktnummer	11055-10
Kategorie	Forceps & Tweezers > Toothed Most Popular > Forceps & Tweezers
EAN	4251909501679



Finde 11055-10 auf finescience.com

Technische Informationen

Versandgewicht	10 g	Nettogewicht	7 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	10 cm		
Tipdimensions	0.8 x 0.7 mm		
Tipwidth	0.8 mm		
Tips	1x2 Teeth		
Tipshape	Curved		