



Graefe Forceps - 1x2 Teeth/Curved/10cm



General Info

Product Number	11055-10
Category	Forceps & Tweezers > Toothed Most Popular > Forceps & Tweezers
EAN	4251909501679



Find 11055-10 on
finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	7 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	10 cm		
Tipdimensions	0.8 x 0.7 mm		
Tipwidth	0.8 mm		
Tips	1x2 Teeth		
Tipshape	Curved		