

F · S · T

## Forceps - Sliding Lock/Smooth/Straight/11cm



## General Info

Product Number	11062-10
Category	Forceps & Tweezers > Smooth
EAN	4251909501686



Find 11062-10 on  
[finescience.com](https://www.finescience.com)

## Technical Info

Shipping Weight	35 g	Net Weight	10 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	11 cm		
Tipdimensions	1 x 0.75 mm		
Tipwidth	1 mm		
Tips	Smooth		
Tipshape	Straight		