

F · S · T

## Delicate Forceps - Smooth/Angled 45°/8cm



The McPherson suture tying forceps are a high-quality precision instrument and used as a ligature aid. The tips form a thread platform with broken edges to prevent the suture thread from being accidentally severed. They are angled at 45° and rounded to minimize damage to the tissue. The spacers in the handle ensure that too much pressure is not built up.

## General Info

Product Number	11063-07
Category	Forceps & Tweezers > Smooth
EAN	4251909501693



Find 11063-07 on  
[finescience.com](https://www.finescience.com)

## Technical Info

Shipping Weight	30 g	Net Weight	6 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	8 cm		
Tipdiameter	0.4 mm		
Tipdimensions	0.4 x 0.5 mm		
Tips	Smooth with platform		
Tipshape	Angled 45°		
Tyingplatform	6 mm		