



Graefe Iris Forceps - 1x2 Teeth/Curved/7cm



General Info

Product Number	11067-07
Category	Forceps & Tweezers > Toothed
EAN	4251909501730



Find 11067-07 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	4 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	7 cm		
Tipdimensions	0.6 x 0.5 mm		
Tipwidth	0.6 mm		
Tips	1x2 Teeth		
Tipshape	Curved		