



Bishop-Harmon Forceps - Serrated/Straight/8.5cm



The fenestrated handle makes these light and ensure a good grip.

General Info

Product Number	11069-08
Category	Forceps & Tweezers > Serrated
EAN	4251909501747



Find 11069-08 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	8.5 cm		
Tipdimensions	0.5 x 0.5 mm		
Tipwidth	0.5 mm		
Tips	Serrated		
Tipshape	Straight		