



Episcleral Forceps - Straight/10cm



The tips of these forceps arch toward each other, closing flat. Exceptional gripping power is possible with minimal damage to tissue. A limit pin prevents excessive applied pressure.

General Info

Product Number	11071-10
Category	Forceps & Tweezers > Smooth
EAN	4251909501761



Find 11071-10 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	10 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	10 cm		
Tipdimensions	1.3 x 3 x 0.2 mm		
Tips	Arched		
Tipshape	Straight		