



Circle Forceps - Smooth/Straight/11.5cm



Cupped forceps for atraumatic gripping of organs and tissues.

General Info

Product Number	11072-11
Category	Forceps & Tweezers > Smooth
EAN	4251909501778



Find 11072-11 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	15 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11.5 cm		
Tipdimensions	7 x 0.5 mm		
Tips	Circle		
Tipshape	Straight		