



Self-Closing Forceps - Cross Serrated/Angled/10cm



These tissue transport forceps with gentle action close flat 4mm wide and 8mm deep.

General Info

|                |                               |
|----------------|-------------------------------|
| Product Number | 11075-00                      |
| Category       | Forceps & Tweezers > Serrated |
| EAN            | 4251909501808                 |



Find 11075-00 on  
[finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |     |
|------------------------|---------------------------------|-----------------|-----|
| Shipping Weight        | 15 g                            | Net Weight      | 9 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA  |
| Minimum Order Quantity | 1                               |                 |     |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |     |
| Length                 | 10 cm                           |                 |     |
| Tipdimensions          | 4 x 8 mm                        |                 |     |
| Tipwidth               | 4 mm                            |                 |     |
| Tips                   | Cross serrated                  |                 |     |
| Tipshape               | Angled                          |                 |     |