




DeBakey Round Handled Forceps



Ideal for gentle, atraumatic gripping of vessels, tissue and other delicate structures. Two rows of atraumatic teeth interfit with one opposite row of atraumatic teeth. The round handles allow easy fingertip adjustments.

General Info

Product Number	11076-15
Category	Forceps & Tweezers > Serrated
EAN	4251909511623



Find 11076-15 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	28 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	15 cm		
Tipwidth	1.2 mm		
Tips	DeBakey serrations		
Tipshape	Straight		