



Flagg Forceps - Serrated/Angled /45°/16cm



General Info

Product Number	11080-02
Category	Forceps & Tweezers > Serrated
EAN	4251909501815



Find 11080-02 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	20 g	Net Weight	16 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	16 cm		
Tipdimensions	1 x 0.9 mm		
Tips	Serrated		
Tipshape	Angled		