



The 1x2 teeth forceps No. 11084-07 are the finest of their kind with 0.1mm wide teeth. The smooth forceps No. 11083-07 have a narrow 5mm long platform. Both forceps have a pin inside the shanks to limit the applied pressure.

General Info

Product Number	11084-07
Category	Forceps & Tweezers > Toothed
EAN	4251909501839



Find 11084-07 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	4 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	7 cm		
Teethwidth	0.1 mm		
Tipdimensions	0.3 x 0.25 mm		
Tipwidth	0.3 mm		
Tips	1x2 Teeth		
Tipshape	Straight		