



Allis Forceps - 5x6 Teeth/Straight/15cm



Tissue and organ holding forceps.

General Info

Product Number	11091-15
Category	Forceps & Tweezers > Toothed
EAN	4251909501860



Find 11091-15 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	29 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	1.5 mm		
Length	15 cm		
Tipdimensions	5 x 3 mm		
Tips	5x6 Teeth		
Tipshape	Straight		