

F · S · T

Castroviejo Suture Tying Forceps - Smooth with platform



Allgemein

Produktnummer	11093-10
Kategorie	Forceps & Tweezers > Smooth
EAN	4251909501884



Finde 11093-10 auf
[finescience.com](https://www.finescience.com)

Technische Informationen

Versandgewicht	35 g	Nettogewicht	12 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	10.5 cm		
Tipdimensions	0.3 x 0.3 mm		
Tipwidth	0.3 mm		
Tips	Smooth with platform		
Tipshape	Straight		
Tyingplatform	5 mm		