



Castroviejo Suture Tying Forceps - Smooth with platform



General Info

|                |                             |
|----------------|-----------------------------|
| Product Number | 11093-10                    |
| Category       | Forceps & Tweezers > Smooth |
| EAN            | 4251909501884               |



Find 11093-10 on  
[finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |      |
|------------------------|---------------------------------|-----------------|------|
| Shipping Weight        | 35 g                            | Net Weight      | 12 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA   |
| Minimum Order Quantity | 1                               |                 |      |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |      |
| Length                 | 10.5 cm                         |                 |      |
| Tipdimensions          | 0.3 x 0.3 mm                    |                 |      |
| Tipwidth               | 0.3 mm                          |                 |      |
| Tips                   | Smooth with platform            |                 |      |
| Tipshape               | Straight                        |                 |      |
| Tyingplatform          | 5 mm                            |                 |      |