



Ring Forceps - with Screw/Smooth/Straight/8.5cm/ID: 1 mm OD: 1.5 mm



For transport of tissue and organs. A limit screw prevents excessive applied pressure.

General Info

Product Number	11101-09
Category	Forceps & Tweezers > Smooth
EAN	4251909501914



Find 11101-09 on finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	10 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	8.5 cm		
Tipdimensions	ID: 1 mm OD: 1.5 mm x 0.5 mm		
Tips	Ring		
Tipshape	Straight		