

F · S · T

Ring Forceps - with Screw/Smooth/Straight/9cm/ID: 2.2 mm OD: 3 mm



For transport of tissue and organs. A limit screw prevents excessive applied pressure.

General Info

| | |
|----------------|-----------------------------|
| Product Number | 11103-09 |
| Category | Forceps & Tweezers > Smooth |
| EAN | 4251909501921 |



Find 11103-09 on
[finescience.com](https://www.finescience.com)

Technical Info

| | | | |
|------------------------|---------------------------------|-----------------|------|
| Shipping Weight | 15 g | Net Weight | 11 g |
| Package Size | 1 Pc. | Unit of Measure | EA |
| Minimum Order Quantity | 1 | | |
| Alloymaterial | Surgical Stainless Steel (Inox) | | |
| Autoclavesafe | Yes | | |
| Length | 9 cm | | |
| Tipdimensions | ID: 2.2 mm OD: 3 mm x 1.3 mm | | |
| Tips | Ring | | |
| Tipshape | Straight | | |