



Dumont #4 Mini Forceps - Standard Tips/Straight/Inox/7cm



General Info

Product Number	11200-12
Category	Forceps & Tweezers > Smooth
EAN	4251909502058



Find 11200-12 on
finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	4 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	7 cm		
Style	Mini #4		
Tipdimensions	0.14 x 0.06 mm		
Tips	Standard		
Tipshape	Straight		