



Dumont #5 Mini Forceps - Standard Tips/Straight/Inox/8.2cm



General Info

Product Number	11200-14
Category	Forceps & Tweezers > Smooth
EAN	4251909502065



Find 11200-14 on
finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	8.2 cm		
Style	Mini #5		
Tipdimensions	0.1 x 0.05 mm		
Tips	Standard		
Tipshape	Straight		