

F · S · T

Dumont SS Forceps - Standard Tips/Straight/Inox/13.5cm



The tips are smooth on the inside for a gentle grip. Applying full strength will not separate tips and a strong hold on the smallest specimen is ensured.

General Info

Product Number	11203-23
Category	Forceps & Tweezers > Smooth
EAN	7640241754732



Find 11203-23 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	15 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloy material	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	13.5 cm		
Style	SS		
Tip dimensions	0.2 x 0.12 mm		
Tips	Standard		
Tip shape	Straight		