



Dumont SS Forceps - Standard Tips/Angled 45°/Inox/13.5cm



The tips are smooth on the inside for a gentle grip. Applying full strength will not separate tips and a strong hold on the smallest specimen is ensured.

General Info

Product Number	11203-25
Category	Forceps & Tweezers > Smooth
EAN	4251909502089



Find 11203-25 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	11 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13.5 cm		
Style	SS		
Tipdimensions	0.2 x 0.12 mm		
Tips	Standard		
Tipshape	Angled 45°		