



Dumont AA Forceps - Standard Tips/Straight/Epoxy Coated/12.5cm



Slip-proof acid and electric shock resistant epoxy coating covers the forceps 17.5mm from the tips.

General Info

Product Number	11210-10
Category	Forceps & Tweezers > Smooth MedTech Line
EAN	7640241754206



Find 11210-10 on finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	17 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Epoxy coated Surgical Stainless Steel (Inox)		
Autoclavesafe	No		
Length	12.5 cm		
Style	AA		
Tipdimensions	0.4 x 0.2 mm		
Tips	Standard		
Tipshape	Straight		