



Dumont AA Forceps - Standard Tips/Straight/Inox/12.5cm



The surface of the tweezers is polished.

General Info

|                |                                             |
|----------------|---------------------------------------------|
| Product Number | 11210-20                                    |
| Category       | Forceps & Tweezers > Smooth<br>MedTech Line |
| EAN            | 4251909502096                               |



Find 11210-20 on  
[finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |      |
|------------------------|---------------------------------|-----------------|------|
| Shipping Weight        | 15 g                            | Net Weight      | 15 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA   |
| Minimum Order Quantity | 1                               |                 |      |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |      |
| Length                 | 12.5 cm                         |                 |      |
| Style                  | AA                              |                 |      |
| Tipdimensions          | 0.4 x 0.2 mm                    |                 |      |
| Tips                   | Standard                        |                 |      |
| Tipshape               | Straight                        |                 |      |