



Dumont #15A Cutter - Carbon Steel/11.5cm



Cuts sutures, fine wire, cartilage, bone, steel wire up to 0.06mm, and copper wire up to 0.15mm. Not autoclavable.

General Info

Product Number	11215-01
Category	Forceps & Tweezers > Cutting
EAN	7640241754305



Find 11215-01 on
finescience.com

Technical Info

Shipping Weight	31 g	Net Weight	30 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Carbon steel		
Autoclavesafe	No		
Bladematerial	Carbon steel		
Cuttingedge	10 mm		
Length	11.5 cm		
Esdsafe	No		
Jawopeningdiameter	2.5 mm		
Tipdimensions	12 x 4.25 mm		
Tips	Cutting		
Tipshape	Angled		