



Dumont #15AGW Cutting Forceps - Epoxy Coated Carbon Steel



High precision cutting tweezers are suitable for applications which require cutting small material or performing small cuts. The strong blades are designed to apply moderate cutting force. Suitable for cut small, soft metal strands (i.e. Copper).

General Info

Product Number	11215-03
Category	Forceps & Tweezers > Cutting MedTech Line
EAN	4251909511890



Find 11215-03 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	26 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Epoxy coated carbon steel		
Autoclavesafe	No		
Bladematerial	Carbon steel		
Cuttingedge	10 mm		
Length	11.5 cm		
Esdsafe	Yes		
Tipdimensions	9 x 4 mm		
Tips	Cutting		
Tipshape	Angled		