



Dumont #15ARW Cutting Forceps - Parallel Blades, Epoxy Coated Carbon Steel



Cuts sutures, fine wire, cartilage, bone, steel wire up to 0.06mm, and copper wire up to 0.15mm. Not autoclavable.

General Info

Product Number	11215-12
Category	Forceps & Tweezers > Cutting MedTech Line
EAN	4251909511906



Find 11215-12 on finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	26 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Epoxy coated carbon steel		
Autoclavesafe	No		
Bladematerial	Carbon steel		
Cuttingedge	7 mm		
Length	11.5 cm		
Esdsafe	Yes		
Jawopeningdiameter	2.5 mm		
Tipdimensions	6 x 4 mm		
Tips	Cutting		
Tipshape	Angled, parallel blades		