



Dumont #2AP Forceps - Standard Tips/Straight/Epoxy Coated/12cm



Provides good insulation for up to 220 AC voltage work and highly impervious to chemical solutions. Frozen or super cold sections will not adhere to the gripping surfaces.

General Info

Product Number	11220-21
Category	Forceps & Tweezers > Smooth MedTech Line
EAN	7640241754251



Find 11220-21 on finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Epoxy coated Surgical Stainless Steel (Inox)		
Autoclavesafe	No		
Length	12 cm		
Style	#2AP		
Tipdimensions	1.5 x 0.7 mm		
Tipwidth	1.5 mm		
Tips	Standard		
Tipshape	Straight		