



Dumont #2 Forceps - Standard Tips/Straight/Inox/12cm



These forceps are ideal for breaking off delicate bones.

General Info

Product Number	11223-20
Category	Forceps & Tweezers > Smooth
EAN	4251909502119



Find 11223-20 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	12 cm		
Style	#2		
Tipdimensions	0.4 x 0.02 mm		
Tips	Standard		
Tipshape	Straight		