



Dumont #2E Components Forceps - 10 x 1.2 mm



General Info

Product Number	11229-10
Category	Forceps & Tweezers > Smooth
EAN	4251909511616



Find 11229-10 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	18 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Dumoxel		
Length	11.5 cm		
Style	#2E		
Tipdimensions	10 x 1.2 mm		
Tips	For components		
Tipshape	Angled		