



Dumont #3C Forceps - Standard Tips/Straight/Inox/11cm



General Info

Product Number	11231-20
Category	Forceps & Tweezers > Smooth MedTech Line
EAN	4251909502126



Find 11231-20 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	11 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11 cm		
Style	#3c		
Tipdimensions	0.17 x 0.1 mm		
Tips	Standard		
Tipshape	Straight		