



Dumont #4 Forceps - Biology Tips/Straight/Dumoxel/11cm



Similar to #5 but with broader shanks.

General Info

Product Number	11242-40
Category	Forceps & Tweezers > Smooth
EAN	4251909502133



Find 11242-40 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Dumoxel		
Antimagnetic	Yes		
Length	11 cm		
Style	#4		
Tipdimensions	0.06 x 0.04 mm		
Tips	Biology		
Tipshape	Straight		