



Dumont #5 Mirror Finish Forceps - Standard Tips/Straight/Inox/11cm



These forceps are highly polished to prevent embryos, nerves, or other adhesive tissues from sticking to the tips.

General Info

Product Number	11251-23
Category	Forceps & Tweezers > Smooth
EAN	4251909502171



Find 11251-23 on finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11 cm		
Style	#5		
Tipdimensions	0.1 x 0.06 mm		
Tips	Standard		
Tipshape	Straight		