



Dumont #5/45 Forceps - Standard Tips/Angled 45°/Dumoxel/11cm



General Info

Product Number	11251-35
Category	Forceps & Tweezers > Smooth MedTech Line Most Popular > Forceps & Tweezers
EAN	7640241754374



Find 11251-35 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Dumoxel		
Antimagnetic	Yes		
Length	11 cm		
Style	#5/45		
Tipdimensions	0.1 x 0.06 mm		
Tips	Standard		
Tipshape	Angled 45°		