



Dumont #5SF Forceps - Super Fine Tips/Straight/Inox/11cm



These forceps have the finest tips ever produced by Dumont.

General Info

Product Number	11252-00
Category	Forceps & Tweezers > Smooth Most Popular > Forceps & Tweezers
EAN	4251909502195



Find 11252-00 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11 cm		
Style	#5SF		
Tipdimensions	0.025 x 0.005 mm		
Tips	Super fine		
Tipshape	Straight		