



Dumont #5 WA Forceps - Biology Tips/Straight/Inox/12cm



These newly designed forceps offer improved stability with an ergonomic body that allows the forceps to fit more comfortably in your hand. A domed end with no soldering grooves allows the forceps to be thoroughly and easily cleaned, preventing biofilm residues.

General Info

Product Number	11252-21
Category	Forceps & Tweezers > Smooth
EAN	4251909502218



Find 11252-21 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	15 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	12 cm		
Style	WA #5		
Tipdimensions	0.05 x 0.02 mm		
Tips	Biology		
Tipshape	Straight		